New Hi-Standard Series

SLIT LAMP SL-95/45DX/45/40

SL-95 SPECIFICATION

MICROSCOPE	Туре	Binocular, Galilean-Type	
	Magnification Changing Method	5 Step Drum Rotation	
	Eyepieces	×12.5 High-Eyepoint	
	Filter for Eyepieces	Yellow Filter	
	Total Magnification	×6.5, ×10, ×16, ×25, ×40	
	Field of View	φ37, φ23,φ14, φ9, φ5.9mm	
	Pupillary Distance	$50 \text{mm} \sim 80 \text{mm}$	
	Diopter Adjustment Range	±7D	
BASE	Working Distance (Back & Forth)	100mm	
	Working Distance (Left & Right)	100mm	
	Working Distance (Up & Down)	30mm	
CHIN REST	Vertial Working Distance	70mm	
	Fixation Target	Tungsten	
ILLUMINATION SYSTEM	Slit Width	0~15mm Continuous	
	Slit Length	0.2, 1, 3, 5, 8, 15mm	
	Light Source (Bulb)	6V 20W (Halogen)	
	Filter	Cobalt, Heat Absorption, Green	
POWER	Voltage	AC 100-120V, 220V, 240V	
	Power Input	50VA	
SIZE	Weight	20~22kg	
	Size	480mm(W) × 410mm(D) × 670mm(H)	

SL-40/45/45DX SPECIFICATION		SL-45DX/45	SL-40	
MICROSCOPE	Туре	Binocular, Galilean-Type	Binocular, Greenough-Type	
	Magnification Changing Method	Drum Rotation 3 Step	Drum Rotation 2 Step	
	Eyepieces	×13.5 High-Eyepoint (×18 available as an option)	×10 High-Eyepoint (×18 available as an option)	
	Total Magnification	×10, ×16, ×25	×10, ×16	
	Field of View	φ20, φ13, φ8mm	φ17, φ10.5, φ6.5mm	
	Pupillary Distance	43mm~80mm	52mm~81mm, 48mm~77mm	
	Diopter Adjustment Range	±5D		
BASE	Working Distance (Back & Forth)	72mm		
	Working Distance (Left & Right)	100mm		
	Working Distance (Up & Down)	30mm		
CHIN REST	Vertial Working Distance	70mm		
	Fixation Target	Tungsten		
ILLUMINATION SYSTEM	Slit Width	0~15mm Continuous		
	Slit Length	0.2, 1, 3, 5, 8, 15mm		
	Light Source (Bulb)	6V 20W (Halogen)		
	Filter	Cobalt, Heat Absorption, Green		
POWER	Voltage	AC 100-120V, 220V, 240V		
	Power Input	50VA		
SIZE	Weight	20~22kg		
	Size	480mm(W) × 410mm(D) × 670mm(H)		

STANDARD ACCESSORIES

■Operation Manual
Power Cable
■Dust Cover
■Chin-rest Paper
■Replacement Halogen Lamp
■Test Bar
■Breath Shield
■Spare Fuse
Replacement Fixation Light Bulb
■Allen Key

OPTIONAL ACCESSORIES

Pro-Pix (Foot Switch Included)*
Tonometer
CCD Camera*
Beam Splitter*
Back Light*
Measuring Eyepiece (x10)**

- * Can not be used for SL-40/45
- **is same manufacturer as the enlisted
- Slit Lamps (None indication is other manufacturer's product)

Manufacturer



TOWA MEDICAL INSTRUMENTS Co., Ltd. 1698-2 Oaza Nakano Nakano-shi, Nagano-ken, 383-0013 Japan

Distributed by



Design and specifications are subject to change without notice.



Contact MEC Sales Division 2-8-4, Kandatsukasa-machi, Chiyoda-ku Tokyo, 101-0048, Japan TEL 81-3-3256-7701 FAX 81-3-3256-7702 E-mail: eye@rexxam.co.jp URL: http://www.rexxam.co.jp URL: http://www.shin-nippon.jp

Printed in Japan I-140101

SHIN-NIPPON by Rexxam

Hi-Standard model to enforce our Slit Lamp line-ups

New Hi-Standard Series SLIT LAMP SL-95 45DX / 45 / 40

New Hi-Standard Series

Complete Slit Lamp line-ups to meet various needs and requirements.

SL-95

In pursuit to user's need of daily examination, SL-95 is a model that improved in optics and illumination system to meet various requirements.

This model is very user friendly and easy to use and also, the digital camera system is attachable for custom usage.

We proudly announce the "New Generation Hi-Standard Model".

Convergence Optics

12degree convergence optics for easy image fusion Newly designed optical alignment affording wider visibility in the evepiece.



Built-In Yellow Filter

Internal yellow filter may be conveniently switched in/out quickly for fluorescein techniques







Reviewed and redesigned the illumination pathway and increased the illumination using the same light bulb. (20W) *Brightness comparison against our existing 20W model.

Expanded Slit Width & Height

A wider slit width & longer slit height capability are attractive features when combined with our new wide field optics Slit width 0-15mm Slit height: 0.2, 1, 3, 5, 8, 15 mm / 0.5 - 15 mm

Diffuser

Standard adjustable diffuser for clear image capturing.





SL-45DX /45 is the standard model that holds every basic function and is also capable to attach digital camera to suit the user's various examination. *Digital camera can only be mounted to SL-45DX



DIGITAL CAMERA SYSTEM

Database Software - Standard with Pro-Pix

Software which can be installed into your PC for Photo and Video of patient examinations. Stored images may be edited



Power Cable Not Required

Pro-Pix Camera is powered by the USB port of your PC, helps to reduce additional cable clutter



Wide Field & High Eye point Eyepiece Design



With our newly created wide field optics, we have successfully increased the actual field of view. Wider fields of view as 37, 23, 14.6, 9.5 & 5.9 mm are achieved. With the new high-point eyepiece feature, eye exams are performed in a more natural, comfortable manner.

Foot Switch



Upon using the digital camera system, the foot switch will let you easily and quickly take pictures of the patient's eye during examination. *Digital Camera System (Pro-Pix DC-200)

is an option.





SL-40 is the basic model that is specialized for every basic function and usability. It is very easy to operate and also capable for quick, precise examination.

Digital Camera System "Pro-Pix DC-200" is available as an option for Model: SL-95 and SL-45DX.